


PERFECT FOR PATIENTS
PERFECT FOR CAREGIVERS



CREATING INNOVATION AND DIFFERENCE

Duet

Duet Mattress System is a medium to high risk system for patients who are at risk of developing Pressure Ulcers according to your sound clinical judgment. Duet can also be used for patients who have an existing stage 1,2,3 and 4 Pressure Ulcer, in conjunction with your hospitals policy on pressure area management.



Various top covers are designed and available to satisfy your individual requirement in terms of therapy, prevention and patient comfort.

Mattress Control Unit

- Dimension: (L)41"x (M)26"x (H)25"
- Weight: 14.9 lbs
- Rated Voltage: 110-120 volts
- Max Current: 5
- Rated Frequency: 60Hz
- Cycle Time: 3~95 mins
- 3 in 1 alternation
- Min/Max Pressure: 15~52 mmHg +/- 4 mmHg
- Static / Pulsation Function
- Protection Type: Class I, Type BF
- 10 Digital Adjustable Comfort Control
- Bottom Out Alarm: Visual / Audio
- Leaving Bed Alarm: Visual / Audio
- Power Failure Alarm: Visual / Audio
- User-Friendly LCD Control Panel
- Fowler Function
- Auto Function
- Turning Function

Mattress Replacement

- Model No: Duet
- Dimension: (L)80"x (W)36"x (H)10"
- Weight: 39.6 lbs
- Weight Capacity: 600lbs
- Cell Material: Blue Nylon with PU Coating
- Cell Number: 18 cells
- Cover Material: Blue Nylon with PU coating finish
- Base Material: Black Polyester with PVC Coating
- CPR for quick CPR release

Toll Free:
877-367-5507

LA, CA OFFICE
TEL: 626-275-4063
TAMPA, FL OFFICE
TEL: 813-854-5464



Duet



Duet


provides you, the carer, with peace of mind
and your patient, with world class care & comfort

Perfect for
patients
and
caregivers

Low slip base cover
is a plus to keep
the mattress
in place
on the bed

The modular
construction
of the cells
makes for simplistic
repair and cleaning

Incorporates with
automatic sensing
technology for
bottoming out
prevention and
pressure
adjustment



PRESSURE MAPPING

Duet provides constant low pressure by redistributing the patient's pressure over a greater area than a standard mattress.


- Equip with 2" convoluted foam base provides additional support during power failure situation

- PU cell design provides easy cleaning and minimizing the proliferation of bacteria

- Mattress replacement can accommodate patient weights up to 600 lbs

- Added safety feature for continuous monitoring patient's status on bottoming out and leaving bed.

- Low Air Loss, Pulsation, 3-1 Alternation, Static and Rotation Therapies in one system



Duet can detect the pressure of patients' buttocks automatically during "Fowler Mode" and adjust the pressure based on these results.

Automatically sets optimal pressure based on patient's weight, height and position.

30 seconds of pulsation simulates a massage to assist in maximizing a patients comfort.

Duet
is your partner
in protection
and treatment
of pressure ulcer

Cable management
keeps
electrical cables
off the floor

CPR quick
release system deflates
the system
within 20 seconds.

Hanging bracket
fits on
all type of beds

The mattress can be adjusted for "Fowler Mode" while acting in conjunction with the bed frame, assisting the medical staff while performing medical procedures.

Multiple lateral rotations provide alternative treatments to patients and medical staffs. Rotation therapy can be set from 3 to 95 minutes.